

BluCapp

CENTRIFUGE TUBES

Designed to cover the everyday need for sample storage and centrifugation in the lab, the **BluCapp Centrifuge tubes** feature easy-to-read graduations as well as a large writing area (20x50mm).

Features:

- The leak-proof caps have large ridges to provide the optimal grip when wearing gloves
- Combined with the 3/4th turn closure, it significantly increases speed and allows one-handed operation
- Made from polypropylene, the **BluCapp Centrifuge tubes** are autoclavable and tested to endure up to 16000 RCF



Ordering Information

Centrifuge Tubes: **BluCapp 15mL and 50mL tubes**

Cat. No.	Description
5100015C	BluCapp 15mL, pre-sterile, 20x25 pcs.
5100050C	BluCapp 50mL, pre-sterile, 20x25 pcs.
5100060C	BluCapp 50mL, pre-sterile, self-standing, 20x25pcs.