


THE ADAPTER TO MOUNT DIGITAL CAMERAS IS PROVIDED TO ALLOW EASY, COST-EFFICIENT DIGITAL IMAGING




**HIGH IMAGES QUALITY AND RICH COLOR**

For various camera types.



**SYSTEM DIAGRAM**



Whether you image for biology, pathology, materials analysis or forensics you will find the adapter is easy to use and provides professional results.

- Easy to setup just only connect the adapter with microscope and camera
- Easy to shoot photographs from microscope
- Easy to control camera via PC, notebook or yourself
- Easy to save images into SD card or your computer directly

The adapter will fit any microscope including compound and stereo microscopes i.e., Olympus, Nikon or Carl Zeiss.

**Relay lens adapter (1001000, 1002000)**



- Use direct to C-Mount on all microscope
- Compatible with DSLR / compact digital camera
  - Canon EOS/ Rebel/ G12, G1X
  - Nikon DSLR/ P7100, P7200

**Eyepiece lens adapter (1007000, 1008000)**



- Use direct eyepiece port on all microscope
- can fit  $\phi$  23.2mm, 30 mm and 30.5 mm eyepiece ports
- Compatible with DSLR camera
  - Canon EOS/ Rebel
- Compatible with Mirrorless DSLR camera
  - Nikon 1 & 2 (model V1, V2, J1, J2 etc)

**HD VDO camera lens adapter (1011000,1012000)**



- Use direct to C-Mount on all microscope
- Compatible with HD VDO camera
  - Canon

**Photo Adapter (1003000,1004000,1005000,1006000)**



- Use direct to photo tube on all microscope
- Compatible with DSLR camera
  - Canon EOS/ Rebel
  - Nikon

**PULSE SCIENCE CO., LTD.**

**Thailand Head Office**

28/9 Moo 1, Bangsrimuang, Muang, Nonthaburi 11000  
 Phone: +66-2886-7808 Fax: 66-886-7809  
 E-mail: [info@pulsescience.co.th](mailto:info@pulsescience.co.th), [sale@pulsescience.co.th](mailto:sale@pulsescience.co.th)

**Singapore Sales Office**

201 Ubi Avenue 4, Singapore 400804  
 Phone: 65-6746-2861 Email: [saleoverseas@pulsescience.co.th](mailto:saleoverseas@pulsescience.co.th)

**Indonesia Sales Office**

**PT Metro Inovasi Asia**  
 Rukan Rose Garden grand Galaxy city blok RRG3 No.31 jakasetia bekasi selatan 17147 Indonesia  
 Phone: 62-81-8888424  
 Email: [sales.indonesia@pulsescience.co.th](mailto:sales.indonesia@pulsescience.co.th), [MetrolnovasiAsia@gmail.com](mailto:MetrolnovasiAsia@gmail.com)

**Website : [www.pulsescience.co.th](http://www.pulsescience.co.th)**