

# IncuWork *ic*



is an IVF Working chamber, totally closed embryo culture system, which substantially enhances the working environment for IVF procedures- ICSI, IMSI, biopsy

The IncuWork ic is designed for efficiency and comfort allowing embryologists to reduce stress and fatigue.

This product can improve embryo implantation rate, clinical pregnancy rate, egg quality and reduce embryo stress.

## FEATURES

- Ergonomic design
- Prepared for inverted microscope for performing ICSI, IMSI, biopsy
- Stable temperature for the entire chamber and hatch box, ambient +5°C to 45°C
- CO2 and O2 level can be controlled with integrated gas mixing system to your preferred level, range 0-10%
- Hatch box: Mechanical interlock is preventing both doors from being opened at the same time.
- Integrated with Anti-Vibration plate
- Touch screen to control and display an overview of all functions
- Alarm system and Data logging
- Remote care concept for check up and maintenance



# IncuWork *iv*



is an IVF work chamber in a closed system. It can help the embryologist reduce the problem of uncontrollable factor of temperature and CO2/O2 during oocytes and embryo handling.

This product provides stable and appropriate environment by precise control of temperature and gas composition which improve the success rate of IVF process.

## FEATURES

- Ergonomic design
- Compatible with all types of stereo microscopes
- Stable temperature throughout the chamber, ambient +5°C to 45°C
- Built-in heating system in the entire working area to maintain a stable temperature of the samples
- CO2 and O2 can be controlled with integrated gas mixing system to your preferred level, range 0-10%
- Easy passing of samples through the side doors and Front Access panel
- Touch screen to control and display an overview of all functions
- Alarm system and Data logging
- Remote care concept for check up and maintenance



## Remote Care Concept



**ADDRESS**  
PULSE SCIENCE CO., LTD.  
28/9 MOO 1 BANGSRIMUANG, MUANG,  
NONTABURI 11000 THAILAND

**TELEPHONE & FAX**  
+66 2886 7808  
+66 2886 7809

**FACEBOOK.COM/PULSESCIENCEGROUP**  
**LINKEDIN.COM/COMPANY/PULSESCIENCE**  
**WWW.PULSESCIENCE.CO.TH**



**ADDRESS**  
PULSE SCIENCE CO., LTD.  
28/9 MOO 1 BANGSRIMUANG, MUANG,  
NONTABURI 11000 THAILAND

**TELEPHONE & FAX**  
+66 2886 7808  
+66 2886 7809

**FACEBOOK.COM/PULSESCIENCEGROUP**  
**LINKEDIN.COM/COMPANY/PULSESCIENCE**  
**WWW.PULSESCIENCE.CO.TH**

