


# IVF LAB MONITORING & ALARM SYSTEM

Solutions for the supply, installation, service, calibration and validation of monitoring & alarm systems to protect your vital products




## Transmitters

Single and multi-channel transmitters capable of measuring and transmitting ;

- Temperature
- Humidity
- Voltage/Current
- Visible Light/UV
- Event
- Differential Pressure
- Gas ; CO<sub>2</sub>/O<sub>2</sub>

## Main Panel Base Station RX250ALD

- 250 channel receiver,
- large memory capacity
- battery backup
- 2 KM Transmissions
- MOXA Communicate with the system over the LAN from multiple PC's



RS232 – USB – Ethernet Connection to PC running **DARCA H1 Software**



- Suite enables users to configure, download and analyses the data transmitted and collected

EMAIL ALARMS — AUTO DIALLER MESSAGE ALARM



ALARM TEXT MESSAGE

- Can text up to 10 phone numbers that can be assigned to specific channels/ items of equipment

## MONITOR UNIT

A standalone sensing unit with easy user recalibration function, remote flashing xenon lamp and sounder, includes local display to show current reading of gas concentration.

- Oxysafe Oxygen Depletion Monitors