

# UDS Scientific

## LABORATORY REFRIGERATORS

---

UDR SERIES



BLOWER SYSTEM (Assisted With Fan)

# UDR-600



Shelf Type



Drawer Type




## GENERAL SPECIFICATIONS

Size		Capacity	Materials			Net weight
External	Internal		External	External door	Internal	
740 x 795 x 1950mm	620 x 635 x 1430mm	610L	Stainless Steel	Glass Door	Stainless Steel	145kg
29.1" x 31.3" x 76.8"	24.4" x 25" x 56.2"	21.5cu.ft				
<b>Cooling system</b>	Single Compressor System	<b>Cooling performance (max)</b>		0°C		
<b>Range of temperature (use)</b>	0°C ~ 10°C	<b>Range of temperature (ambient)</b>		15°C ~ 35°C		
<b>Temperature sensor</b>	PT100 Sensor	<b>Refrigerant</b>		R134a		
<b>Compressor</b>	Hermetic Type	<b>Type of condenser</b>		Block Condenser (Pin, Filter)		
<b>Controller</b>	Digital	<b>Shelf materials &amp; Q'ty</b>	<b>Shelf Type</b>	Included 4 Wire Shelves		
			<b>Drawer Type</b>	Included 7 Stainless Drawers		
<b>Rated voltage</b>	220V, 50Hz, 2.7A	<b>Power consumption</b>		398W		

\* Voltage and Hz can be customized according to customer's needs

\* All specifications are subject to change prior notice

## DETAIL VIEW



BLOWER SYSTEM (Assisted With Fan)

# UDR-1125-1C



Shelf Type





Drawer Type



CUF - SUF & CCF SERIES

FD & SCF SERIES

OPTIONS & ACCESSORIES

UD SERIES

UDR SERIES

BLOWER SYSTEM (Assisted With Fan)

# UDR-1125-2C



Shelf Type





Drawer Type



CUF - SUF & CCF SERIES

FD & SCF SERIES

OPTIONS & ACCESSORIES

UD SERIES

UDR SERIES

## GENERAL SPECIFICATIONS

Size		Capacity	Materials			Net weight	
External	Internal		External	Internal door	Internal	UDR-1125-1C	UDR-1125-2C
1260 x 795 x 1950mm	1140 x 635 x 1430mm	1125L	Stainless Steel	Glass Door	Stainless Steel	207kg	230kg
49.6" x 31.3" x 76.8"	44.8" x 25" x 56.2"	39.7cu.ft					
Cooling system	UDR-1125-1C	Single Compressor System	Cooling performance (max)			0°C	
	UDR-1125-2C	Dual compressor system	Range of temperature (ambient)			15°C - 35°C	
Range of temperature (use)		0°C - 10°C	Refrigerant			R134a	
Temperature sensor		PT100 Sensor	Type of condenser			Block Condenser (Pin, Filter)	
Compressor		Hermetic Type	Shelf materials & Q'ty	Shelf Type		Included 8 Wire Shelves	
Controller		Digital		Drawer Type		Included 12 Stainless Drawers	
Rated voltage	UDR-1125-1C	220V, 50Hz, 3.2A	Power consumption	UDR-1125-1C		536W	
	UDR-1125-2C	220V, 50Hz, 5.8A		UDR-1125-2C		980W	

\* Voltage and Hz can be customized according to customer's needs

\* All specifications are subject to change prior notice

## DETAIL VIEW

