



Semen Quality Analyzer

Compact & Automatic



1 Easy & fast

LensHooke™ is easy to use. No technical background is needed to perform the test. Results are obtained in 2-5 minutes.



2 Touch Panel

LensHooke™ uses a resistive touch screen to allow medical personnel to operate screen with gloves on.



3 Auto Focus & Full HD

LensHooke™ is capable of taking 1080P Full HD dynamic images to enable immediate sperm morphology examination by medical personnel.



4 Real-Time Image

LensHooke™ is equipped with an HDMI connection for external monitors. It can display dynamic images with all key analytical data (pH/concentration /motility/morphology) on screen in real time.

4 Key Data for Diagnosis



pH value



Concentration



Morphology



Motility
(NP, PR, Total Motility)

※ Comply with WHO 4th & WHO 5th edition.

Operation Steps

